



Recombinant HEV Non-structural polyprotein (genotype 4) (DAGCF-0165)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Recombinant 6xHis tagged HEV Non-structural polyprotein (genotype 4) (UniProtKB No. Q91VZ9) was expressed.
Purity	>85%
Conjugate	His
Molecular Weight	187Kda
Recommended Usage	This item requires custom production and lead time is between 4-114 weeks. We can custom produce based on your specifications.
Format	lyophilized
Size	1 mg, 5 mg
Buffer	Each vial contains 100 µg of purified protein in PBS
Preservative	None
Storage	Store at 2-8 °C short term and -80°C long term. It is recommended to aliquot needed volumes and store at -80°C. Avoid repeated freeze/thaw cycles.

BACKGROUND

Keywords HCV;HCV Non-structural;Hepatitis C Virus Non-structural;Non-structural;Hepatitis C virus

GENE INFORMATION

