



# Recombinant HPV16 E6 protein [His] (DAGF-094)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Antigen Description</b>	Recombinant Papillomavirus 16 E6 antigen produced in E.Coli contains 159 amino covering the full length of HPV-16 having its C-terminus fused with a 6xHis tag and a total molecular weight of 21kDa. The HPV16 was purified by proprietary chromatographic technique.
<b>Purity</b>	>95%
<b>Conjugate</b>	His
<b>Molecular Weight</b>	21kDa
<b>Format</b>	Liquid
<b>Concentration</b>	Batch dependent - please inquire should you have specific requirements.
<b>Size</b>	100 µg
<b>Buffer</b>	25mM Tris-Base, 25mM k2CO3 and 1M Urea
<b>Preservative</b>	None
<b>Storage</b>	Recombinant HPV-16 E6 although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.

## BACKGROUND

<b>Keywords</b>	HPV 16; HPV; Human papillomavirus (HPV); Papillomaviridae; Human Papilloma Virus Type 16; Human Papilloma Virus; Human papillomavirus type 16 E6; Human papillomavirus type 16 E6
-----------------	---