



Pyrethroid [BSA] (DAGF-203)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Pyrethroids are axonic excitotoxins, the toxic effects of which are mediated through preventing the closure of the voltage-gated sodium channels in the axonal membranes. The sodium channel is a membrane protein with a hydrophilic interior.
Nature	Synthetic
Expression System	Synthetic
Conjugate	BSA
Procedure	None
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	0.01M pH7.4 PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Keywords Methothrin; MB pyrethroid; Methothryn; 34388-29-9; 4-(Methoxymethyl)benzyl chrysanthemum monocarboxylate; 4-(Methoxymethyl)benzyl chrysanthemum monocarboxylate