



HIV type 1 gp120 (IIIB, Lyophilized) [FITC] (DAGH002)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Baculovirus derived gp120 HIV type 1 IIIB conjugated to FITC
Specificity	This fluoresce conjugated gp120 binds to murine monoclonal antibodies of defined epitope specificity
Conjugate	FITC
Applications	Dot Blots; FACS; Immunohistocytology; Direct Immunofluorescence
Stability	6 months at -80°C
Bio-activity	This fluoresce conjugated gp120 binds to CD4 bearing T-cells as determined by FACS.
Format	Lyophilized. Solubilize in DI H2O for use. Store solubilized gp120 at -20°C for long term storage
Size	50 T/50 µg
Preservative	None
Storage	Room Temp

BACKGROUND

Keywords	HIV type 1 gp120; Human Immunodeficiency Virus-1 gp120; HIV-1 gp120; Envelope glycoprotein GP120; gp120
-----------------	---