



# Recombinant HIV type 1 P24 Protein [FITC] (DAGH004)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Antigen Description	Recombinant p24 HIV type 1 covalently coupled to FITC of 0.7nm diameter.
Specificity	FITC-p24 binds to Murine mAbs of defined epitope specificity and Human pAbs as determined by dot blot assays.
Conjugate	FITC
Applications	T-Cell Research; Diagnostics
Stability	At least 12 Months at 4°C
Bio-activity	ND
Size	50 µg
Preservative	None
Storage	-75 °C

## BACKGROUND

Keywords	HIV type 1 p24; HIV-1 p24; HIV1 p24; CA; Capsid protein p24; Human immunodeficiency virus 1; Human immunodeficiency virus type 1 p24; Retroviridae; Lentivirus
----------	--