



## HIV type 1 vif (IIIB) [FITC] (DAGH008)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

Antigen Description	E.coli derived vif (HIV type 1 ) conjugated to FITC.
Specificity	This fluoresce conjugated vif binds to murine monoclonal antibodies of defined epitope specificity
Conjugate	FITC
Applications	Dot Blots; FACS; Immunohistocytology; Direct Immunofluorescence
Stability	At least 6 months at -75°C
Bio-activity	Not determined
Format	Frozen Liquid
Size	50 µg
Preservative	None
Storage	-75 °C

### BACKGROUND

Keywords	HIV type 1 vif; Viral infectivity factor; HIV-1 vif; HIV1 vif; Human immunodeficiency virus 1; Human immunodeficiency virus type 1 vif; Human immunodeficiency virus 1 Viral infectivity factor
----------	---