



HIV type 1 nef (IIIB) [FITC] (DAGH009)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	E.coli derived nef (HIV type 1 IIIB) conjugated to FITC.
Specificity	This fluoresce conjugated nef binds to murine monoclonal antibodies of defined epitope specificity.
Conjugate	FITC
Applications	Dot Blots; FACS; Immunohistocytology; Direct Immunofluorescence
Stability	At least 6 months at -75°C
Bio-activity	Not determined
Format	Frozen Liquid
Size	50 µg
Preservative	None
Storage	-75 °C

BACKGROUND

Keywords	HIV type 1 nef; HIV-1 Nef; Human immunodeficiency virus type 1 Nef; Negative Regulatory Factor; Protein Nef; 3'ORF; Negative factor
-----------------	---