



HIV type 1 gp120 (YU2) [FITC] (DAGH011)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Baculovirus derived gp120 HIV type 1 YU2 (M-tropic) conjugated to FITC.
Specificity	This fluoresce conjugated gp120 binds to murine monoclonal antibodies of defined epitope specificity.
Conjugate	FITC
Applications	Dot Blots; FACS; Immunohistocytology; Direct Immunofluorescence
Stability	6 months at -75°C
Bio-activity	This fluoresce conjugated gp120 binds to CD4 bearing T-cells as determined by FACS
Format	Frozen Liquid
Size	50 µg
Preservative	None
Storage	-75 °C

BACKGROUND

Keywords	HIV type 1 gp120; Human Immunodeficiency Virus-1 gp120; HIV-1 gp120; Envelope glycoprotein GP120; gp120
-----------------	---