



HIV type 1 gp120 (ADA) [FITC] (DAGH013)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Baculovirus derived gp120 (HIV type 1 ADA) conjugated to FITC.
Specificity	This fluoresce conjugated gp120 binds to murine monoclonal antibodies of defined epitope specificity and immobilized CD4 in Dot Blot assays.
Nature	Recombinant
Expression System	Insect cells
Conjugate	FITC
Applications	Dot Blots; FACS; Immunohistocytology; Direct Immunofluorescence
Stability	At least 6 months at -75°C
Bio-activity	This fluoresce conjugated gp120 binds to CD4 bearing cells in FACS.
Procedure	None
Format	Frozen Liquid
Size	50 µg
Preservative	None
Storage	-75 °C

BACKGROUND

Keywords	HIV type 1 gp120; Human Immunodeficiency Virus-1 gp120; HIV-1 gp120; Envelope glycoprotein GP120; gp120
-----------------	---