



Estriol-3 [HRP] (DAGP-2524H)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Estriol(E3), also spelled oestriol, is a steroid, a weak estrogen, and a minor female sex hormone. It is one of three major endogenous estrogens, the others being estradiol and estrone. Levels of estriol in women who are not pregnant are almost undetectable.
Species	N/A
Conjugate	HRP
Applications	ELISA
Format	Liquid, Linking Position-3-CME position.
Size	1 ml
Buffer	Supplied in Tris with Enzyme Stabilizers.
Preservative	Anti-Bacterial Agent.
Storage	Store at 2-8 °C.