



Recombinant Human T-Cell Lymphotropic Virus Type I gp21 antigen (DAGRM071T)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

| | |
|--------------|---|
| Conjugate | Unconjugated |
| Applications | Single and chimeric recombinant antigen for application as coating in indirect method. For CILAs & TRFIA assays |
| Preservative | None |
| Storage | Store at +4°C for short term (1-2 weeks) and -20°C for long term. Avoid multiple freeze/thaw cycles. |

BACKGROUND

| | |
|--------------|--|
| Introduction | Human T-lymphotropic virus (HTLV) belongs to the family Retroviridae, genus Deltaretrovirus. It is a human, single-stranded RNA retrovirus that causes T-cell leukemia and T-cell lymphoma in adults and may also be involved in certain demyelinating diseases. Adult T-lymphotropic virus (ATLV) is a strain of this disease that effect primarily adults. |
| Keywords | HTLV-1; gp21; gp46 |