



Sterigmatocystin standard (DAGS010)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Use this reference material as calibrant for chromatography or any other analytical technique.
Conjugate	N/A
Cas.No	10048-13-2
Molecular Weight	324.28g/mol
Size	5 mg
Storage	The reference material should be stored in the original sealed container at the $-20\pm 4^{\circ}\text{C}$. After use the bottle should be tightly closed and protected from moisture.
Ship	Ship at room temperature.