



## Nicarbazine standard (DAGS084)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

Antigen Description	Use this reference material as calibrant for chromatography or any other analytical technique.
Expression System	N/A
Purity	98.79%(Please refer to the certificate of analysis of each batch.)
Conjugate	N/A
Cas.No	330-95-0
Molecular Weight	426.38g/mol
Format	Neat
Size	100 mg
Storage	The reference material should be stored in the original sealed container at the 20±4°C. After use the bottle should be tightly closed and protected from moisture.
Ship	Ship at room temperature.