



Magic™ Anti-Rotavirus Monoclonal antibody, Clone C208N (DCAB-TJ118)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Reacts with intact virus of strains RRV, WA, SA-11 and bovine. Immunostains infected cells. Does not react with Influenza A, Influenza B, RSV, Parainfluenza 1, 2 & 3, Adenovirus, M. pneumonia, H. pylori and Mammalian cells.
Target	Rotavirus
Immunogen	Purified Rhesus Rotavirus (RRV)
Isotype	IgG2a
Source/Host	Mouse
Species Reactivity	Rotavirus
Clone	C208N
Purification	> 90% pure. Protein A Chromatography
Conjugate	Unconjugated
Applications	Suitable for use in ELISA and IFA. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded. Recommended pairs fo Suggested pair for testing (Capture - Detection): DCAB-TJ117 - DCAB-TJ118
Format	Purified, Liquid
Concentration	100 µg/mL (OD280 nm, E0.1% = 1.3)
Size	1 mg

Buffer	0.01 M PBS, pH 7.2 This product contains no stabilizing protein.
Preservative	0.1% Sodium Azide
Storage	Upon receipt, aliquot and store at -20°C. Avoid multiple freeze/thaw cycles.

BACKGROUND

Introduction	Rotavirus is the most common cause of severe diarrhea among infants and young children. It is a genus of double-stranded RNA virus in the family Reoviridae. Nearly every child in the world has been infected with rotavirus at least once by the age of five. Immunity develops with each infection, so subsequent infections are less severe; adults are rarely affected. There are five species of this virus, referred to as A, B, C, D, and E. Rotavirus A, the most common species, causes more than 90% of rotavirus infections in humans.
Keywords	Rotavirus