



Anti-AGGF1 monoclonal antibody (DCABH-10455)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	The AGGF1 gene encodes a potent angiogenic factor that contains a forkhead-associated domain and a G-patch domain (Tian et al., 2004 [PubMed 14961121]).
Immunogen	A synthetic peptide of human AGGF1 is used for rabbit immunization.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Protein A
Conjugate	Unconjugated
Applications	Western Blot (Transfected lysate); ELISA
Buffer	In 1x PBS, pH 7.4
Preservative	None
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	AGGF1 angiogenic factor with G patch and FHA domains 1 [Homo sapiens]
Official Symbol	AGGF1
Synonyms	AGGF1; angiogenic factor with G patch and FHA domains 1; FLJ10283; GPATC7; GPATCH7;

HSU84971; VG5Q; angiogenic factor VG5Q; vasculogenesis gene on 5q protein; G patch domain-containing protein 7; HUS84971;

Entrez Gene ID	55109
Protein Refseq	NP_060516
UniProt ID	Q8N302
Chromosome Location	5q13.3
Function	eukaryotic cell surface binding; nucleic acid binding; protein binding;