



Anti-APBB2 monoclonal antibody (DCABH-10541)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	The protein encoded by this gene interacts with the cytoplasmic domains of amyloid beta (A4) precursor protein and amyloid beta (A4) precursor-like protein 2. This protein contains two phosphotyrosine binding (PTB) domains, which are thought to function in signal transduction.
Immunogen	A synthetic peptide of human APBB2 is used for rabbit immunization.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Protein A
Conjugate	Unconjugated
Applications	Western Blot (Transfected lysate); ELISA
Buffer	In 1x PBS, pH 7.4
Preservative	None
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	APBB2 amyloid beta (A4) precursor protein-binding, family B, member 2 [Homo sapiens]
Official Symbol	APBB2

Synonyms	APBB2; amyloid beta (A4) precursor protein-binding, family B, member 2; amyloid beta A4 precursor protein-binding family B member 2; Fe65 like; FE65L; FE65L1; MGC35575; Fe65-like 1; protein Fe65-like 1; DKFZp434E033;
Entrez Gene ID	323
Protein Refseq	NP_001159522
UniProt ID	B4DJ88
Chromosome Location	4p13
Function	beta-amyloid binding; protein binding; transcription factor binding;