



Anti-ARL9 monoclonal antibody (DCABH-10612)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	ARL9 is a member of the small GTPase protein family with a high degree of similarity to ARF (MIM 103180) proteins of the RAS superfamily.
Immunogen	A synthetic peptide of human ARL9 is used for rabbit immunization.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Protein A
Conjugate	Unconjugated
Applications	Western Blot (Transfected lysate); ELISA
Buffer	In 1x PBS, pH 7.4
Preservative	None
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	ARL9 ADP-ribosylation factor-like 9 [Homo sapiens]
Official Symbol	ARL9
Synonyms	ARL9; ADP-ribosylation factor-like 9;

Entrez Gene ID [387029](#)

Chromosome Location 4q12
