



Anti-C4ORF6 monoclonal antibody (DCABH-10791)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	This gene is expressed in neuroblastoma; however, the function of this gene is not yet determined.
Immunogen	A synthetic peptide of human C4orf6 is used for rabbit immunization.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Protein A
Conjugate	Unconjugated
Applications	Western Blot (Transfected lysate); ELISA
Buffer	In 1x PBS, pH 7.4
Preservative	None
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	C4orf6 chromosome 4 open reading frame 6 [Homo sapiens]
------------------	---

Official Symbol	C4orf6
Synonyms	C4ORF6; chromosome 4 open reading frame 6; uncharacterized protein C4orf6; aC1; expressed in neuroblastoma;
Entrez Gene ID	10141
Protein Refseq	NP_005741
Chromosome Location	4p16