



Anti-CINP monoclonal antibody (DCABH-11022)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	The function of this gene has not been reported. This gene contains a trinucleotide (ctt) repeat length polymorphism; alleles with either one or two repeat units have been identified.
Immunogen	A synthetic peptide of human CINP is used for rabbit immunization.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Protein A
Conjugate	Unconjugated
Applications	Western Blot (Transfected lysate); ELISA
Buffer	In 1x PBS, pH 7.4
Preservative	None
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	CINP cyclin-dependent kinase 2 interacting protein [Homo sapiens]
Official Symbol	CINP
Synonyms	CINP; cyclin-dependent kinase 2 interacting protein; cyclin-dependent kinase 2-interacting

protein; MGC849; CDK2-interacting protein;

Entrez Gene ID	51550
Protein Refseq	NP_116019
UniProt ID	A0A024R6M9
Chromosome Location	14q32.33
Function	protein binding;