



Anti-CKAP2 monoclonal antibody (DCABH-11026)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	CKAP2 is a cytoskeleton-associated protein involved in mitotic progression (Seki and Fang, 2007 [PubMed 17376772]).
Immunogen	A synthetic peptide of human CKAP2 is used for rabbit immunization.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Protein A
Conjugate	Unconjugated
Applications	Western Blot (Transfected lysate); ELISA
Buffer	In 1x PBS, pH 7.4
Preservative	None
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	CKAP2 cytoskeleton associated protein 2 [Homo sapiens]
Official Symbol	CKAP2
Synonyms	CKAP2; cytoskeleton associated protein 2; cytoskeleton-associated protein 2; FLJ10749; LB1;

se20 10; TMAP; CTCL tumor antigen se20-10; tumor- and microtubule-associated protein;
tumor-associated microtubule-associated protein; se20-10; DKFZp686L1238;

Entrez Gene ID	26586
----------------	-----------------------

Protein Refseq	NP_001091995
----------------	------------------------------

UniProt ID	Q8WWK9
------------	------------------------

Chromosome Location	13q14
---------------------	-------
