



Anti-CPNE4 monoclonal antibody (DCABH-11117)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Calcium-dependent membrane-binding proteins may regulate molecular events at the interface of the cell membrane and cytoplasm. This gene is one of several genes that encode a calcium-dependent protein containing two N-terminal type II C2 domains and an integrin A domain-like sequence in the C-terminus.
Immunogen	A synthetic peptide of human CPNE4 is used for rabbit immunization.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Protein A
Conjugate	Unconjugated
Applications	Western Blot (Transfected lysate); ELISA
Buffer	In 1x PBS, pH 7.4
Preservative	None
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	CPNE4 copine IV [Homo sapiens]
Official Symbol	CPNE4

Synonyms	CPNE4; copine IV; copine 8; COPN4; CPN4;
Entrez Gene ID	8902
Chromosome Location	3q22.1