



Anti-DIRAS2 monoclonal antibody (DCABH-11286)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	DIRAS2 belongs to a distinct branch of the functionally diverse Ras (see HRAS; MIM 190020) superfamily of monomeric GTPases.
Immunogen	A synthetic peptide of human DIRAS2 is used for rabbit immunization.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Protein A
Conjugate	Unconjugated
Applications	Western Blot (Transfected lysate); ELISA
Buffer	In 1x PBS, pH 7.4
Preservative	None
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	DIRAS2 DIRAS family, GTP-binding RAS-like 2 [Homo sapiens]
Official Symbol	DIRAS2
Synonyms	DIRAS2; DIRAS family, GTP-binding RAS-like 2; GTP-binding protein Di-Ras2; Di Ras2;

DKFZp761C07121; distinct subgroup of the Ras family member 2; Di-Ras2;

Entrez Gene ID	54769
Protein Refseq	NP_060064
UniProt ID	Q96HU8
Chromosome Location	9q22.32
Function	GTP binding; GTPase activity; nucleotide binding;