



Anti-DNAJC12 monoclonal antibody (DCABH-11301)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	This gene encodes a member of a subclass of the HSP40/DnaJ protein family. Members of this family of proteins are associated with complex assembly, protein folding, and export. Two transcript variants encoding distinct isoforms have been identified for this gene.
Immunogen	A synthetic peptide of human DNAJC12 is used for rabbit immunization.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Protein A
Conjugate	Unconjugated
Applications	Western Blot (Transfected lysate); ELISA
Buffer	In 1x PBS, pH 7.4
Preservative	None
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	DNAJC12 DnaJ (Hsp40) homolog, subfamily C, member 12 [Homo sapiens]
Official Symbol	DNAJC12

Synonyms	DNAJC12; DnaJ (Hsp40) homolog, subfamily C, member 12; dnaJ homolog subfamily C member 12; J domain protein 1; JDP1; j domain-containing protein 1; J domain containing protein 1 (JDP1); RP11-57G10.2;
Entrez Gene ID	56521
Protein Refseq	NP_068572
UniProt ID	Q6IAH1
Chromosome Location	10q21.3
Function	heat shock protein binding; unfolded protein binding;