



Anti-DOM3Z monoclonal antibody (DCABH-11316)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	This gene localizes to the major histocompatibility complex (MHC) class III region on chromosome 6. The function of its protein product is unknown, but its ubiquitous expression and conservation in both simple and complex eukaryotes suggests that this may be a housekeeping gene.
Immunogen	A synthetic peptide of human DOM3Z is used for rabbit immunization.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Protein A
Conjugate	Unconjugated
Applications	Western Blot (Transfected lysate); ELISA
Buffer	In 1x PBS, pH 7.4
Preservative	None
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	DOM3Z dom-3 homolog Z (C. elegans) [Homo sapiens]
Official Symbol	DOM3Z

Synonyms	DOM3Z; dom-3 homolog Z (C. elegans); DOM 3 (C. elegans) homolog Z; protein Dom3Z; NG6; RAI1; DOM3L;
Entrez Gene ID	1797
Protein Refseq	NP_005501
UniProt ID	A0A024RCW8
Chromosome Location	6p21.3
Function	metal ion binding; nucleotide binding;