



Anti-E2F7 monoclonal antibody (DCABH-11352)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	E2F transcription factors, such as E2F7, play an essential role in the regulation of cell cycle progression (Di Stefano et al., 2003 [PubMed 14633988]).
Immunogen	A synthetic peptide of human E2F7 is used for rabbit immunization.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Protein A
Conjugate	Unconjugated
Applications	Western Blot (Transfected lysate); ELISA
Buffer	In 1x PBS, pH 7.4
Preservative	None
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	E2F7 E2F transcription factor 7 [Homo sapiens]
Official Symbol	E2F7
Synonyms	E2F7; E2F transcription factor 7; transcription factor E2F7; E2F-7; FLJ12981;

Entrez Gene ID	144455
Protein Refseq	NP_976328
UniProt ID	Q96AV8
Chromosome Location	12q21.1
Pathway	E2F transcription factor network, organism-specific biosystem;
Function	DNA binding; identical protein binding; sequence-specific DNA binding transcription factor activity;