



Anti-EPB41L4A monoclonal antibody (DCABH-11416)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Members of the band 4.1 protein superfamily, including EPB41L4A, are thought to regulate the interaction between the cytoskeleton and plasma membrane (Ishiguro et al., 2000 [PubMed 10874211]).
Immunogen	A synthetic peptide of human EPB41L4A is used for rabbit immunization.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Protein A
Conjugate	Unconjugated
Applications	Western Blot (Transfected lysate); ELISA
Buffer	In 1x PBS, pH 7.4
Preservative	None
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	EPB41L4A erythrocyte membrane protein band 4.1 like 4A [Homo sapiens]
Official Symbol	EPB41L4A

Synonyms	EPB41L4A; erythrocyte membrane protein band 4.1 like 4A; band 4.1-like protein 4A; NBL4; erythrocyte protein band 4.1-like 4; EPB41L4; FLJ38738; DKFZp566L203;
Entrez Gene ID	64097
Protein Refseq	NP_071423
UniProt ID	Q8N8X1
Chromosome Location	5q21.3
Function	binding; cytoskeletal protein binding; molecular_function;