



Anti-EPS8L3 monoclonal antibody (DCABH-11439)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	This gene encodes a protein that is related to epidermal growth factor receptor pathway substrate 8 (EPS8), a substrate for the epidermal growth factor receptor. The function of this protein is unknown. Alternatively spliced transcript variants encoding different isoforms exist.
Immunogen	A synthetic peptide of human EPS8L3 is used for rabbit immunization.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Protein A
Conjugate	Unconjugated
Applications	Western Blot (Transfected lysate); ELISA
Buffer	In 1x PBS, pH 7.4
Preservative	None
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	EPS8L3 EPS8-like 3 [Homo sapiens]
Official Symbol	EPS8L3

Synonyms	EPS8L3; EPS8-like 3; epidermal growth factor receptor kinase substrate 8-like protein 3; FLJ21522; MGC16817; EPS8-like protein 3; EPS8-related protein 3; epidermal growth factor receptor pathway substrate 8-related protein 3; EPS8R3;
Entrez Gene ID	79574
Protein Refseq	NP_078802
UniProt ID	A8K2J6
Chromosome Location	1p13.2