



Anti-FAF2 monoclonal antibody (DCABH-11487)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	The protein encoded by this gene is highly expressed in peripheral blood of patients with atopic dermatitis (AD), compared to normal individuals. It may play a role in regulating the resistance to apoptosis that is observed in T cells and eosinophils of AD patients.
Immunogen	A synthetic peptide of human FAF2 is used for rabbit immunization.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Protein A
Conjugate	Unconjugated
Applications	Western Blot (Transfected lysate); ELISA
Buffer	In 1x PBS, pH 7.4
Preservative	None
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	FAF2 Fas associated factor family member 2 [Homo sapiens]
Official Symbol	FAF2

Synonyms	FAF2; Fas associated factor family member 2; UBX domain containing 8 , UBXD8; FAS-associated factor 2; ETEA; expressed in T cells and eosinophils in atopic dermatitis; KIAA0887; UBX domain protein 3B; UBXN3B; UBX domain containing 8; UBX domain-containing protein 8; UBX domain-containing protein 3B; expressed in T-cells and eosinophils in atopic dermatitis; UBXD8;
Entrez Gene ID	23197
Protein Refseq	NP_055428
UniProt ID	Q96CS3
Chromosome Location	5q35.2
Function	protein binding; ubiquitin binding; ubiquitin protein ligase binding;