



Anti-FBRS monoclonal antibody (DCABH-11509)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Fibrosin is a lymphokine secreted by activated lymphocytes that induces fibroblast proliferation (Prakash and Robbins, 1998 [PubMed 9809749]).
Immunogen	A synthetic peptide of human FBRS is used for rabbit immunization.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Protein A
Conjugate	Unconjugated
Applications	Western Blot (Transfected lysate); ELISA
Buffer	In 1x PBS, pH 7.4
Preservative	None
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	FBRs fibrosin [Homo sapiens]
Official Symbol	FBRs
Synonyms	FBRs; fibrosin; FBS1, fibrosin 1; probable fibrosin-1; FBS; FLJ11618; fibrosin 1; fibrogenic

lymphokine; probable fibrosin-1 long transcript protein; FBS1;

Entrez Gene ID	64319
----------------	-----------------------

Protein Refseq	NP_001098549
----------------	------------------------------

UniProt ID	Q9HAH7
------------	------------------------

Chromosome Location	16p11.2
---------------------	---------
