



Anti-FGFBP2 monoclonal antibody (DCABH-11565)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	This gene encodes a member of the fibroblast growth factor binding protein family. The encoded protein is a serum protein that is selectively secreted by cytotoxic lymphocytes and may be involved in cytotoxic lymphocyte-mediated immunity. An increase in the amount of gene product may be associated with atopic asthma and mild extrinsic asthma. Staff
Immunogen	A synthetic peptide of human FGFBP2 is used for rabbit immunization.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Protein A
Conjugate	Unconjugated
Applications	Western Blot (Transfected lysate); ELISA
Buffer	In 1x PBS, pH 7.4
Preservative	None
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	FGFBP2 fibroblast growth factor binding protein 2 [Homo sapiens]
Official Symbol	FGFBP2
Synonyms	FGFBP2; fibroblast growth factor binding protein 2; fibroblast growth factor-binding protein 2; killer specific secretory protein of 37 kDa; KSP37; FGF-BP2; FGFBP-2; HBp17-RP; FGF-binding protein 2; HBp17-related protein; 37 kDa killer-specific secretory protein; killer-specific secretory protein of 37 kDa; HBP17RP;
Entrez Gene ID	83888
Protein Refseq	NP_114156
UniProt ID	Q9BYJ0
Chromosome Location	4p15.32
Function	growth factor binding;