



Anti-FKBPL monoclonal antibody (DCABH-11579)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	The protein encoded by this gene has similarity to the immunophilin protein family, which play a role in immunoregulation and basic cellular processes involving protein folding and trafficking. The encoded protein is thought to have a potential role in the induced radioresistance. Also it appears to have some involvement in the control of the cell cycle.
Immunogen	A synthetic peptide of human FKBPL is used for rabbit immunization.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Protein A
Conjugate	Unconjugated
Applications	Western Blot (Transfected lysate); ELISA
Buffer	In 1x PBS, pH 7.4
Preservative	None
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	FKBPL FK506 binding protein like [Homo sapiens]
Official Symbol	FKBPL

Synonyms	FKBPL; FK506 binding protein like; FK506-binding protein-like; DIR1; NG7; FK506-binding protein like; WAF-1/CIP1 stabilizing protein 39; WISP39;
Entrez Gene ID	63943
Protein Refseq	NP_071393
UniProt ID	A0A024RCQ6
Chromosome Location	6p21.3
Function	FK506 binding; binding; peptidyl-prolyl cis-trans isomerase activity;