



Anti-FTSJ3 monoclonal antibody (DCABH-11625)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Although the function of this gene is not known, the existence of this gene is supported by mRNA and EST data. A possible function of the encoded protein can be inferred from amino acid sequence similarity to the E.coli FtsJ protein and to a mouse protein possibly involved in embryogenesis.
Immunogen	A synthetic peptide of human FTSJ3 is used for rabbit immunization.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Protein A
Conjugate	Unconjugated
Applications	Western Blot (Transfected lysate); ELISA
Buffer	In 1x PBS, pH 7.4
Preservative	None
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	FTSJ3 FtsJ homolog 3 (E. coli) [Homo sapiens]
Official Symbol	FTSJ3
Synonyms	FTSJ3; FtsJ homolog 3 (E. coli); putative rRNA methyltransferase 3; protein ftsJ homolog 3; rRNA (uridine-2-O-)-methyltransferase 3; EPCS3; FLJ20062;
Entrez Gene ID	117246
Protein Refseq	NP_060117
UniProt ID	Q8IY81
Chromosome Location	17q23
Function	methyltransferase activity; nucleic acid binding; transferase activity;