



Anti-NAP1L3 monoclonal antibody (DCABH-12509)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	This gene encodes a member of the nucleosome assembly protein (NAP) family. It is linked closely to a region of genes responsible for several X-linked mental retardation syndromes. Alternatively spliced transcript variants have been found for this gene, but the biological validity of the variants has not been determined.
Immunogen	A synthetic peptide of human NAP1L3 is used for rabbit immunization.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Protein A
Conjugate	Unconjugated
Applications	Western Blot (Transfected lysate); ELISA
Buffer	In 1x PBS, pH 7.4
Preservative	None
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	NAP1L3 nucleosome assembly protein 1-like 3 [Homo sapiens]
Official Symbol	NAP1L3

Synonyms	NAP1L3; nucleosome assembly protein 1-like 3; MB20; MGC26312; NPL3;
Entrez Gene ID	4675
Protein Refseq	NP_004529
UniProt ID	Q8IYV1
Chromosome Location	Xq21.3-q22