



Anti-APOOL monoclonal antibody (DCABH-14617)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	This gene encodes a protein which contains an apolipoprotein O superfamily domain. This domain is found on proteins in circulating lipoprotein complexes.
Immunogen	A synthetic peptide of human APOOL is used for rabbit immunization.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Protein A
Conjugate	Unconjugated
Applications	Western Blot (Transfected lysate); ELISA
Buffer	In 1x PBS, pH 7.4
Preservative	None
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	APOOL apolipoprotein O-like [Homo sapiens]
Official Symbol	APOOL
Synonyms	APOOL; apolipoprotein O-like; chromosome X open reading frame 33 , CXorf33, FAM121A,

family with sequence similarity 121A; AAIR8193; UNQ8193; family with sequence similarity 121A; CXorf33; FAM121A; MGC129748;

Entrez Gene ID [139322](#)

Protein Refseq [NP_940852](#)

UniProt ID [Q6UXV4](#)

Chromosome Location Xq21.1
