



Anti-CCDC12 monoclonal antibody, clone 8C2 (DCABH-14888)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	CCDC12 (NP_653317.1, 1 a.a. ~ 166 a.a) full-length recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Isotype	IgG2a
Source/Host	Mouse
Species Reactivity	Human
Clone	8C2
Conjugate	Unconjugated
Applications	Western Blot (Cell lysate); Western Blot (Recombinant protein); Sandwich ELISA (Recombinant protein); ELISA
Sequence Similarities	MEATTAGVGRLEEEALRRKERLKALEKTGRKDKEDGEPKTKHLREEEEEEGEKHRELRLR NYVPEDEDLKKRRVPQAKPVAVEEKVKEQLEAAKPEPVIEEVDLANLAPRKPDWDLKRDV AKKLEKLKKRTQRAIAELIRERLKGQEDSLASAVDAATEQKTCDS
Buffer	In 1x PBS, pH 7.4
Preservative	None
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name [CCDC12 coiled-coil domain containing 12 \[Homo sapiens \]](#)

Official Symbol	CCDC12
Synonyms	CCDC12; coiled-coil domain containing 12; coiled-coil domain-containing protein 12; MGC23918; FLJ39430; FLJ40801;
Entrez Gene ID	151903
Protein Refseq	NP_653317
UniProt ID	Q8WUD4
Chromosome Location	3p21.31
Pathway	Spliceosome, organism-specific biosystem; Spliceosome, conserved biosystem; Spliceosome, 35S U5-snRNP, organism-specific biosystem.