



# Anti-CCDC87 monoclonal antibody (DCABH-14906)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Immunogen</b>	A synthetic peptide of human CCDC87 is used for rabbit immunization.
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Purification</b>	Protein A
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	Western Blot (Transfected lysate); ELISA
<b>Buffer</b>	In 1x PBS, pH 7.4
<b>Preservative</b>	None
<b>Storage</b>	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

## GENE INFORMATION

<b>Gene Name</b>	<a href="#">CCDC87 coiled-coil domain containing 87 [ Homo sapiens ]</a>
<b>Official Symbol</b>	CCDC87
<b>Synonyms</b>	CCDC87; coiled-coil domain containing 87; coiled-coil domain-containing protein 87; FLJ10786;
<b>Entrez Gene ID</b>	<a href="#">55231</a>

<b>Protein Refseq</b>	<a href="#">NP_060689</a>
<b>UniProt ID</b>	<a href="#">Q9NVE4</a>
<b>Chromosome Location</b>	11q13.2