



Anti-LRGUK monoclonal antibody, clone 8C5 (DCABH-15955)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	LRGUK (NP_653249, 721 a.a. ~ 825 a.a) partial recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Isotype	IgG2a
Source/Host	Mouse
Species Reactivity	Human
Clone	8C5
Conjugate	Unconjugated
Applications	Western Blot (Cell lysate); Western Blot (Recombinant protein); ELISA
Sequence Similarities	YFKPPFGPYPEKSGKDSLVSMLKCSLFRFCPWSKELPFQPPEGSISSHLGSGASDSETEET RKALPIQSFSHEKESHQHRQHSVPVISRPGSNVKPTLPPIQGR
Buffer	In 1x PBS, pH 7.4
Preservative	None
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	LRGUK leucine-rich repeats and guanylate kinase domain containing [Homo sapiens]
Official Symbol	LRGUK

Synonyms	LRGUK; leucine-rich repeats and guanylate kinase domain containing; leucine-rich repeat and guanylate kinase domain-containing protein; FLJ32786;
Entrez Gene ID	136332
Protein Refseq	NP_653249
UniProt ID	Q96M69
Chromosome Location	7q33
Function	ATP binding; kinase activity; nucleotide binding; transferase activity;