



Anti-TBC1D22B monoclonal antibody (DCABH-17050)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	TBC1D22B may act as a GTPase activating protein for Rab family protein(s).
Immunogen	A synthetic peptide of human TBC1D22B is used for rabbit immunization.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Protein A
Conjugate	Unconjugated
Applications	Western Blot (Transfected lysate); ELISA
Buffer	In 1x PBS, pH 7.4
Preservative	None
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	TBC1D22B TBC1 domain family, member 22B [Homo sapiens]
Official Symbol	TBC1D22B
Synonyms	TBC1D22B; TBC1 domain family, member 22B; C6orf197, chromosome 6 open reading frame 197; TBC1 domain family member 22B; dJ744I24.2; FLJ20337; C6orf197; RP4-744I24.2;

MGC125626; MGC125627; DKFZp762J0110;

Entrez Gene ID	55633
-----------------------	-----------------------

Protein Refseq	NP_060242
-----------------------	---------------------------

UniProt ID	Q9NU19
-------------------	------------------------

Chromosome Location	6p21.2
----------------------------	--------

Function	GTPase activator activity; Rab GTPase activator activity;
-----------------	---
