



Anti-PBX3 monoclonal antibody, clone FQS3944(3) (DCABH-176)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit monoclonal to Pbx3
Antigen Description	Transcriptional activator that binds the sequence 5-ATCAATCAA-3.
Immunogen	Synthetic peptide (the amino acid sequence is considered to be commercially sensitive)
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse, Rat, Human
Clone	FQS3944(3)
Purity	Tissue culture supernatant
Conjugate	Unconjugated
Applications	WB
Positive Control	Human fetal brain, A431 and HeLa lysates.
Format	Liquid
Size	100 µl
Buffer	PBS 49%,Sodium azide 0.01%,Glycerol 50%,BSA 0.05%
Preservative	0.1% Sodium Azide
Storage	Store at -20°C. Stable for 12 months at -20°C

GENE INFORMATION

Gene Name	PBX3 pre-B-cell leukemia homeobox 3 [Homo sapiens]
Official Symbol	PBX3
Synonyms	PBX3; pre-B-cell leukemia homeobox 3; pre B cell leukemia transcription factor 3; pre-B-cell leukemia transcription factor 3; homeobox protein PBX3;
Entrez Gene ID	5090
Protein Refseq	NP_001128250
UniProt ID	P40426
Chromosome Location	9q33.3
Pathway	Transcriptional misregulation in cancer, organism-specific biosystem; Transcriptional misregulation in cancer, conserved biosystem;
Function	sequence-specific DNA binding; sequence-specific DNA binding transcription factor activity;