



Anti-ZNF446 monoclonal antibody (DCABH-17614)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	May be involved in transcriptional regulation.
Immunogen	A synthetic peptide of human ZNF446 is used for rabbit immunization.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Protein A
Conjugate	Unconjugated
Applications	Western Blot (Transfected lysate); ELISA
Buffer	In 1x PBS, pH 7.4
Preservative	None
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	ZNF446 zinc finger protein 446 [Homo sapiens]
Official Symbol	ZNF446
Synonyms	ZNF446; zinc finger protein 446; FLJ20626; ZKSCAN20; zinc finger protein with KRAB and SCAN domains 20; ZSCAN30;

Entrez Gene ID	55663
Protein Refseq	NP_060378
UniProt ID	Q9NWS9
Chromosome Location	19q13.43
Pathway	Gene Expression, organism-specific biosystem; Generic Transcription Pathway, organism-specific biosystem.
Function	DNA binding; metal ion binding; sequence-specific DNA binding transcription factor activity; zinc ion binding;