



Anti-GAS2L2 monoclonal antibody (DCABH-200611)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	The protein encoded by this gene appears to crosslink microtubules and microfilaments and may be part of the cytoskeleton. This gene is mainly expressed in skeletal muscle.
Immunogen	A synthetic peptide of human GAS2L2 is used for rabbit immunization.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Protein A
Conjugate	Unconjugated
Applications	WB, ELISA
Size	1 mg
Buffer	In 1x PBS, pH 7.4
Preservative	None
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	GAS2L2 growth arrest-specific 2 like 2 [Homo sapiens (human)]
Official Symbol	GAS2L2
Synonyms	GAS2L2; growth arrest-specific 2 like 2; GAR17; GAS2-like protein 2; GAS2-related protein on chromosome 17;
Entrez Gene ID	246176
Protein Refseq	NP_644814
UniProt ID	Q8NHY3
Chromosome Location	17q12
Function	actin filament binding; cytoskeletal adaptor activity; microtubule binding;