



Anti-CCSER1 monoclonal antibody (DCABH-200958)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	CCSER1 (coiled-coil serine-rich protein 1) is a protein-coding gene. Diseases associated with CCSER1 include schizophrenia.
Immunogen	A synthetic peptide of human MGC48628 is used for rabbit immunization.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Protein A
Conjugate	Unconjugated
Applications	WB, ELISA
Size	1 mg
Buffer	In 1x PBS, pH 7.4
Preservative	None
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	CCSER1 coiled-coil serine-rich protein 1 [Homo sapiens (human)]
Official Symbol	CCSER1

Synonyms	CCSER1; coiled-coil serine-rich protein 1; FAM190A; serine-rich coiled-coil domain-containing protein 1; protein FAM190A; family with sequence similarity 190, member A;
Entrez Gene ID	401145
Protein Refseq	NP_001138537
UniProt ID	Q9C0I3
Chromosome Location	4q22.1