



Anti-NOL7 monoclonal antibody (DCABH-201052)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	NOL7 (nucleolar protein 7, 27kDa) is a protein-coding gene. Diseases associated with NOL7 include cervical cancer, and cervicitis.
Immunogen	A synthetic peptide of human NOL7 is used for rabbit immunization.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Protein A
Conjugate	Unconjugated
Applications	WB, ELISA
Size	1 mg
Buffer	In 1x PBS, pH 7.4
Preservative	None
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	NOL7 nucleolar protein 7, 27kDa [Homo sapiens (human)]
Official Symbol	NOL7

Synonyms	NOL7; nucleolar protein 7, 27kDa; PQBP3; RARG-1; C6orf90; dJ223E5.2; nucleolar protein 7; nucleolar protein of 27 kDa; polyglutamine binding protein 3; retinoic acid repressible protein;
Entrez Gene ID	51406
Protein Refseq	NP_057251
UniProt ID	A0A024QZW2
Chromosome Location	6p23
Function	poly(A) RNA binding;