



Anti-PLEKHN1 monoclonal antibody (DCABH-201270)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	PLEKHN1 (pleckstrin homology domain containing, family N member 1) is a protein-coding gene. GO annotations related to this gene include phospholipid binding.
Immunogen	A synthetic peptide of human PLEKHN1 is used for rabbit immunization.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Protein A
Conjugate	Unconjugated
Applications	WB, ELISA
Size	1 mg
Buffer	In 1x PBS, pH 7.4
Preservative	None
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	PLEKHN1 pleckstrin homology domain containing, family N member 1 [Homo sapiens (human)]
------------------	---

Official Symbol	PLEKHN1
Synonyms	PLEKHN1; pleckstrin homology domain containing, family N member 1; pleckstrin homology domain-containing family N member 1; PH domain-containing family N member 1;
Entrez Gene ID	84069
Protein Refseq	NP_001153656
UniProt ID	Q494U1
Chromosome Location	1p36.33