



Anti-THAP2 monoclonal antibody (DCABH-201578)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	THAP2 (THAP domain containing, apoptosis associated protein 2) is a protein-coding gene. GO annotations related to this gene include DNA binding and metal ion binding. An important paralog of this gene is THAP3.
Immunogen	A synthetic peptide of human THAP2 is used for rabbit immunization.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Protein A
Conjugate	Unconjugated
Applications	WB, ELISA
Size	1 mg
Buffer	In 1x PBS, pH 7.4
Preservative	None
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name [THAP2 THAP domain containing, apoptosis associated protein 2 \[Homo sapiens \(human\) \]](#)

Official Symbol	THAP2
Synonyms	THAP2; THAP domain containing, apoptosis associated protein 2; THAP domain-containing protein 2;
Entrez Gene ID	83591
Protein Refseq	NP_113623
UniProt ID	Q9H0W7
Chromosome Location	12q21.1
Function	DNA binding; metal ion binding;