



Anti-SLAMF6 monoclonal antibody, clone FQS0027 (DCABH-4571)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit monoclonal to LY108
Antigen Description	Triggers cytolytic activity only in natural killer cells (NK) expressing high surface densities of natural cytotoxicity receptors.
Immunogen	Synthetic peptide corresponding to residues within Human LY108 (UniProt Q96DU3).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Clone	FQS0027
Conjugate	Unconjugated
Applications	WB, IP
Positive Control	Ramos; Jurkat and MOLT4 cell lysates.
Format	Liquid
Size	100 µl
Buffer	Preservative: 0.01% Sodium azide; Constituents: 50% Glycerol, 0.05% BSA
Preservative	0.01% Sodium Azide
Storage	Store at -20°C.

Ship

Shipped at 4°C.

GENE INFORMATION

Gene Name	SLAMF6 SLAM family member 6 [Homo sapiens]
Official Symbol	Slamf6
Synonyms	SLAMF6; SLAM family member 6; CD352; KALI; KAL1b; Ly108; NTB A; NTBA; SF2000; NTBA receptor; NK-T-B-antigen; activating NK receptor; natural killer-, T- and B-cell antigen; NTB-A; FLJ50657; MGC104953;
Entrez Gene ID	114836
Protein Refseq	NP_001171643
UniProt ID	Q96DU3
Chromosome Location	1q23.1
Function	receptor activity;