



Anti-SSSCA1 monoclonal antibody, clone 3G6 (DCABH-564)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Mouse monoclonal to SSSCA1
Antigen Description	Might play a role in mitosis. Antigenic molecule. Could be a centromere-associated protein. May induce anti-centromere antibodies.
Immunogen	Recombinant full length Human SSSCA1 produced in HEK293T cells (NP_006387).
Isotype	IgG1
Source/Host	Mouse
Species Reactivity	Human
Clone	3G6
Purification	This antibody is purified from Mouse ascites fluid by affinity chromatography.
Conjugate	Unconjugated
Applications	WB, Flow Cyt, ICC/IF
Positive Control	HEK293T cells and cell lysate transfected with pCMV6-ENTRY SSSCA1 cDNA; HeLa and Jurkat cells; COS7 cells transiently transfected by pCMV6-ENTRY SSSCA1.
Format	Liquid
Size	100 µl
Buffer	pH: 7.30; Preservative: 0.02% Sodium azide; Constituents: 48% PBS, 1% BSA, 50% Glycerol
Preservative	0.02% Sodium Azide

Storage	store at -20°C. Avoid repeated freeze / thaw cycles.
----------------	--

Ship	Shipped at 4°C.
-------------	-----------------

GENE INFORMATION

Gene Name	SSSCA1 Sjogren syndrome/scleroderma autoantigen 1 [Homo sapiens]
------------------	--

Official Symbol	SSSCA1
------------------------	--------

Synonyms	SSSCA1; Sjogren syndrome/scleroderma autoantigen 1; Sjogrens syndrome/scleroderma autoantigen 1; Sjogren syndrome/scleroderma autoantigen 1; p27; autoantigen p27; centromeric autoantigen (27kD);
-----------------	--

Entrez Gene ID	10534
-----------------------	-----------------------

Protein Refseq	NP_006387
-----------------------	---------------------------

UniProt ID	O60232
-------------------	------------------------

Chromosome Location	11q13.1
----------------------------	---------

Function	identical protein binding;
-----------------	----------------------------