



Anti-TTLL12 monoclonal antibody, clone 2E4 (DCABH-656)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Mouse monoclonal to TTLL12
Immunogen	Recombinant full length Human TTLL12 produced in HEK293T cells (NP_055955).
Isotype	IgG1
Source/Host	Mouse
Species Reactivity	Human
Clone	2E4
Purification	This antibody was purified from mouse ascites fluids by affinity chromatography.
Conjugate	Unconjugated
Applications	WB, Flow Cyt
Positive Control	HEK293T cells transfected with the pCMV6-ENTRY TTLL12
Format	Liquid
Size	100 µl
Buffer	pH: 7.30; Preservative: 0.02% Sodium azide; Constituents: 48% PBS, 1% BSA, 50% Glycerol
Preservative	0.02% Sodium Azide
Storage	store at -20°C. Avoid freeze / thaw cycles.
Ship	Shipped at 4°C.

GENE INFORMATION

Gene Name	TTLL12 tubulin tyrosine ligase-like family, member 12 [Homo sapiens]
Official Symbol	TTLL12
Synonyms	TTLL12; tubulin tyrosine ligase-like family, member 12; tubulin--tyrosine ligase-like protein 12; KIAA0153; dJ526I14.2; FLJ41795;
Entrez Gene ID	23170
Protein Refseq	NP_055955
UniProt ID	A0A024R4U3
Chromosome Location	22q13.31
Function	tubulin-tyrosine ligase activity;