



# Anti-FIBP monoclonal antibody (DCABH-687)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

|                            |  |
|----------------------------|--|
| <b>Product Overview</b>    | Mouse monoclonal to FIBP   |
| <b>Antigen Description</b> | May be involved in mitogenic function of FGF1.   |
| <b>Immunogen</b>           | Recombinant full length protein, corresponding to amino acids 1-358 of Human FIBP with Proprietary tag (tag alone is 26 KDa) (AAH14388). |
| <b>Isotype</b>             | IgG2a  |
| <b>Source/Host</b>         | Mouse  |
| <b>Species Reactivity</b>  | Human  |
| <b>Conjugate</b>           | Unconjugated   |
| <b>Applications</b>        | Sandwich ELISA   |
| <b>Positive Control</b>    | Recombinant FIBP protein   |
| <b>Format</b>              | Liquid   |
| <b>Size</b>                | 100 µg   |
| <b>Buffer</b>              | pH: 7.20; Constituent: 99% PBS   |
| <b>Preservative</b>        | None   |
| <b>Storage</b>             | Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.   |

## GENE INFORMATION

**Gene Name** [FIBP fibroblast growth factor \(acidic\) intracellular binding protein \[ Homo sapiens \]](#)

|                            |  |
|----------------------------|--|
| <b>Official Symbol</b>     | FIBP   |
| <b>Synonyms</b>            | FIBP; fibroblast growth factor (acidic) intracellular binding protein; acidic fibroblast growth factor intracellular-binding protein; FGFIBP; aFGF intracellular-binding protein; FGF-1 intracellular-binding protein; FIBP-1; |
| <b>Entrez Gene ID</b>      | <a href="#">9158</a>   |
| <b>Protein Refseq</b>      | <a href="#">NP_004205</a>  |
| <b>UniProt ID</b>          | <a href="#">O43427</a>   |
| <b>Chromosome Location</b> | 11q13.1  |
| <b>Function</b>            | fibroblast growth factor binding;  |