



Anti-CCDC115 monoclonal antibody (DCABH-689)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Mouse monoclonal to MGC12981
Immunogen	Recombinant fragment, corresponding to amino acids 49-180 of Human MGC12981 (AAH06429.1).
Isotype	IgG2b
Source/Host	Mouse
Species Reactivity	Human
Conjugate	Unconjugated
Applications	WB, IHC-P, Sandwich ELISA
Positive Control	Recombinant MGC12981 protein; MGC12981 transfected 293T cell lysate; Human breast tissue.
Format	Liquid
Size	100 µg
Buffer	pH: 7.20; Constituent: 99% PBS
Preservative	None
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	CCDC115 coiled-coil domain containing 115 [Homo sapiens]
Official Symbol	CCDC115
Synonyms	CCDC115; coiled-coil domain containing 115; coiled-coil domain-containing protein 115; ccp1; FLJ30131; MGC12981;
Entrez Gene ID	84317
Protein Refseq	NP_115733
UniProt ID	A0A024QZZ6
Chromosome Location	2q21.1